

Stability and Shipping information for iScreen Liquid Control Urine

STORAGE & STABILITY CLAIMS– (Reprinted from package insert)

Unopened:

- A. The controls are stable until the expiration date when stored at -10 to -20°C and protected from light.
- B. The controls are stable until the expiration date when stored at 2-8°C, however, no stability claims can be made for Oxazepam as it may deteriorate over time when stored refrigerated.

After Opening:

- A. The controls are stable for six months or until the expiration date, whichever comes first, when stored at -10 to -20°C. (Controls can be aliquoted and frozen)
- B. The controls are stable for 31 days or until the expiration date, whichever comes first, when stored tightly capped at 2-8°C.
- C. Thaw controls as needed; allow to come to room temperature followed by gentle swirling before use.

LIMITATIONS CONCERNING THC AND OXAZEPAM:

Controls should be stored unopened either frozen or refrigerated until expiration date. Long term storage is best frozen, as the Oxazepam value can drop when stored refrigerated for longer than 10 months. THC will bind to certain plastics and should not be aliquoted into plastic micro tubes or regular plastic test tubes.

ISCREEN CONTROLS, OXAZEPAM STABILITY: (Reprinted from package insert)

Oxazepam has known stability problems in urine stored refrigerated, and to a lesser degree, frozen. Our experience indicates that Oxazepam will not deteriorate more than 10% of target level for at least one year when stored frozen at -20°C. Further deterioration may occur beyond this period although Oxazepam ordinarily tests positive throughout the control's shelf life.

ISCREEN CONTROLS, THC STABILITY: (Abstracted from package insert)

iScreen controls are stable for the length of time under the storage conditions stated in the package insert. In spite of this fact, under certain conditions there may be observed a gradual decline in THC levels from repeated sampling of a single bottle of control material. This drop in THC values may occur from any THC sample (i.e. calibrators, controls, and samples). This decline of THC most often occurs from handling and not from product instability. It is well known that THC-COOH binds to surfaces, especially certain plastics.

ALIUQUOTING CONTROLS:

If end-user desires the control can safely undergo 5 thaw/freeze cycles without affecting performance. The control can be thawed and returned to freezer to extend shelf life to 6 months. Thaw controls as needed by allowing to come to room temperature followed by gentle swirling before use.

SHIPPING METHODS:

If shipped within the continental United States, iScreen controls are frozen then shipped one day with ice packs.

If shipped internationally, Frozen iScreen controls are shipped with ice packs in insulated coolers.

The controls will not remain frozen upon customer receipt. End-user should promptly place controls in either freezer or refrigerator. The controls do not need to remain frozen, but excessive heat (>35°C) will decrease Oxazepam 10% or greater depending on length of time heated.