

Strep A Test Kit

Preview Strep A Test Kit is a rapid chromatographic immunoassay intended for use by healthcare professionals as a qualitative *in vitro* diagnostic test for detection of group A streptococcal antigen from throat swab specimens.




- Easy-to-use
- Convenient all-in-one test kit
- CLIA-waived
- Room temperature storage
- Results in approximately 10 minutes
- Test and treat during the same office visit




Test Procedure

1 Specimen Collection




2 Step 1


Add 4 drops Add 4 drops



3 Step 2




4 Step 3




Immerse test strip to the level of the marker line.

5 Step 4

5-10 MIN



6 Results



Positive Negative Invalid

Control line → C Test line → T

Ordering information

Item Code	Description	Kit Size / Case Size	Complexity
PRE-STA-25	25-Strip; Strep A Throat Swab Test	25 / 20	CLIA WAIVED